


PROTOCOL & PRODUCT OVERVIEW



MTP RT

RNA-seq workflow integration of Multi-Targeted Priming (MTP) Reverse Transcription (RT)

-  Seamless protocol integration
-  No specialized equipment necessary
-  Combined target selection and off-target counter-selection
-  Genome-wide scalability



The MTP RT
technology



Organism-Specific Solution



We use annotation data to
design custom solutions for
each organism.

All we need is a RefSeq ID.

Flexibly
Defined
Targeting

Target Sequence Space

30,000+ unique targets

Typically:

✓ mRNA

Can flexibly include:

✓ Non-polyadenylated mRNA

✓ lncRNA

✓ miRNA

✓ Pathogen transcripts

Off-Target Sequence Space

5,000+ unique off-targets

Typically:

✗ rRNA

✗ tRNA

Can flexibly include:

✗ Globin mRNA

✗ Uninformative mRNAs


✗ Host background

Intelligent
Oligo
Design

Graph Representation Learning
Enrichment Platform

Parallelized
Oligo
Synthesis

Proprietary Parallelized
Oligo Synthesis



Sample Preparation

Reverse Transcription

Library Preparation

Bulk & sc RNA-seq



Seamless
Protocol
Integration



Multi-Targeted Priming
Reagent



Compared to Oligo(dT), MTP gets...

Low prevalence
transcripts

31–331% more

31-331% increased number of well
measured genes

Differential
expression

66–136% more

66-136% increased number of identified
significantly differentially expressed
genes (GEL₅₀ threshold)

Genome-wide
enrichment

30,000+ targets

Unbiased enrichment across gene pathways

No more 3' bias

Non-PolyA transcripts

Lose the 3' bias and target hard-to-reach
mitochondrial RNA, miRNA, lncRNA, viral
RNA and prokaryotic mRNA, and any other
non-polyadenylated transcripts

Statistics compiled from bulk RNA sequencing, microarray profiling, and qPCR
results for *S. cerevisiae* and *N. crassa*.



MTP Applications

Non-Model Organisms

Use our online product catalog or search by exact RefSeq Assembly ID at multitarga.com/search to find the non-model organism MTP offerings. These MTPs target the mRNA transcripts while counter-selecting for the rRNA and tRNA using the available sequences in RefSeq.

Model Organisms

Our model-organism MTPs target mRNAs while counter-selecting for rRNAs and tRNAs, using their RefSeq and/or GenBank assemblies as the sequence source material. Contact us for additional customization options including globin counter-selection, target gene-set restriction, and more.



Human MTP



Arabidopsis MTP



Mouse MTP



Drosophila MTP




Rat MTP



Zebrafish MTP

Host-Pathogen

Adaptively define target and off-target sequence spaces to enable counter-selection of background host material. Specific host-pathogen offerings coming soon. Contact us for custom collaboration on your host-pathogen or other multi-organism samples.



Pan-Onco

Comprehensive enrichment of 769 cancer-associated genes across all tissue types for precision oncology transcriptomics. Check out our website, multitarga.com/pan-onco, to find the full target gene list.

769

Cancer genes

Microbiome

We are currently developing reverse transcription and capture solutions for the human microbiome and other microbial communities. Contact us to collaborate on our metagenomic and metatranscriptomic solutions.



Overview

Multi-Targeted Priming (MTP) is a novel reverse transcription strategy that uses flexibly designed primer pools to enrich transcripts of interest. This approach augments the first step of bulk RNA-seq library preparation, reverse transcription, by replacing traditional oligo(dT) or random hexamer primers with our context specific reagent.

MTP is fully compatible with existing laboratory equipment and reagents. It is integrated seamlessly into the reverse transcription step for any Bulk RNA-seq library preparation kit without requiring changes to enzymes, buffers, or thermal cycler programs.

No new instrumentation required. MTP is a one-step swap-in replacement for randomers or oligo(dT) primers during cDNA synthesis.

Key Benefits

- Increased transcript specificity during cDNA synthesis
- Drop-in compatibility with existing RNA-seq protocols and kits
- Effective with total RNA, FFPE, and degraded samples
- No changes to equipment or enzymes
- **Cost-effective enrichment** prior to library prep
- **Flexible input amounts** (10–500 ng total RNA)
- **Supplied as lyophilized oligo pools** for long-term storage



Compatible Library Preparation Kits

MTPs can be used with any RNA-seq library preparation method that includes a first-strand cDNA synthesis step. Below are validated or compatible:

Kit Name	Company	Cat. No.	Integration Description
SuperScript IV First-Strand Synthesis System	Thermo Fisher	18091050	Targeted primers allowed for RT
Maxima H Minus First Strand cDNA Synthesis Kit	Thermo Fisher	K1651	Targeted primers allowed for RT
ProtoScript II First Strand cDNA Synthesis Kit	NEB	E6560	Designed for targeted primer usage
TruSeq Stranded mRNA	Illumina	20020594	Supports upstream RT prior to library prep
TruSeq Stranded Total RNA with Ribo-Zero	Illumina	20020596	Supports upstream RT prior to library prep
RNA Prep with Enrichment	Illumina	20040536	Compatible with user-prepared cDNA workflows
NEBNext Ultra II RNA Library Prep	NEB	E7770	Accepts externally Prepared cDNA
NEBNext Ultra II Directional RNA Library Prep	NEB	E7760	Accepts externally Prepared cDNA
AmpliSeq for Illumina Custom RNA Panel	Illumina	20020496	RT step separable; Randomer substitution for MTPs

Contact us if your kit is not listed. We are validating protocol integration.



Protocol: MTP Reverse Transcription

Input

- RNA input: 10–500 ng of high-quality or partially degraded total RNA
- DNase treatment to remove residual genomic DNA
- Compatible with FFPE RNA and rRNA-depleted samples

Step 1: Annealing of MTP Primers

1. In a nuclease-free PCR tube, combine the following components:

Component	Volume	Notes
Total RNA	Up to 11 μ L	Adjust to 11 μ L total with water
MTP primer (2 μ M)	1 μ L	Replace oligo(dT) or N6 primers
dNTP mix (10 mM each)	1 μ L	Final 0.5 mM concentration per nucleotide

2. Mix gently by pipetting or flicking the tube
3. Spin briefly to collect liquid at the bottom
4. Place in a thermal cycler and heat to $T_m + 2^\circ\text{C}$ for 5 minutes (refer to the MTP vial label for T_m)
5. Immediately place the tube on ice for at least 1 minute

Step 2: First-Strand Reverse Transcription Reaction

1. Add the following to the RNA/primer mix (final volume 20 μ L):

Component	Volume	Notes
5X First-Strand RT buffer	4 μ L	Use buffer supplied with your RT enzyme
100 mM DTT	1 μ L	Reducing agent for optimal enzyme activity
Rnase Inhibitor (40 U/ μ L)	1 μ L	e.g., RNaseOUT or equivalent
Reverse Transcriptase	1 μ L	e.g., SuperScript IV, ProtoScript II, etc.

2. Mix gently and briefly spin down
3. Incubate in a thermal cycler using the following program
 - 42–50°C for 50 minutes (based on enzyme recommendation)
 - 70°C for 15 minutes to inactivate the enzyme
 - Hold at 4°C or proceed immediately to downstream steps



MTP Reverse Transcription - Protocol

Continued.

Step 3: Post-RT Workflow

- Follow the library prep workflow with the resulting cDNA.
- Continue with second-strand synthesis, adapter ligation, tagmentation, depletion, or enrichment as required by your selected kit.
- MTP-generated cDNA is fully compatible with stranded workflows, dUTP-based marking, and downstream hybrid capture.

MTP Storage & Stability

Form: Supplied as lyophilized (dried) oligo pool reagent

Storage: Store at -20°C in a desiccated container

Resuspension: Use nuclease-free water or 10 mM Tris-HCl pH 8.0

Aliquoting: Recommended to avoid freeze-thaw

Shelf life: Stable >1 year at -20°C

Reaction Yield

MTP vial quantity	Total Oligo Yield	Estimated RT Reactions (20 μL , 2 μM MTP)
25 nmole	~25,000 pmol	250-500 reactions
100 nmole	~100,000 pmol	1,000-2,000 reactions

Contact

For questions, assistance, or comments email us at contact@multitarga.com.
<https://multitarga.com>.